





DreamStar™ Intro

DreamStar™Intro Evolve w/ water reservoir

# DreamStar™ Intro DreamStar™ Intro Evolve

DreamStar™ Intro is designed to improve the patient's compliance with the therapy by providing performance, comfort and ease of use. This makes it easily adaptable to the needs of the user.

#### **COMPACT DESIGN**

DreamStar™ Evolve as a CPAP device has also a fully integrated heated humidifier increasing the respiratory comfort, while remaining one of the most compact devices in the market. The humidifier is located at the rear of the device with its water reservoir. The reservoir and the air inlet filter are easy to remove for cleaning.

The fully integrated power supply and heated humidifier, makes DreamStar™ Evolve one of the lightest devices in the market. This ensures maximum comfort while using the device either traveling or at home.

### **CONFIGURATIONS**

The DreamStar™ Intro has a configuration without humidifier and the Evolve configuration with an integrated heated humidifier. The Evolve Configuration may be used either with a water reservoir for humidification or without humidification by replacing the reservoir with a cover.

## Advanced DreamStar™ Water & Humid Control Technology

The innovative design of the water reservoir "DreamStar™ Water Control" and the automatic humidification system "Humid Control" provide a new concept for water management. This prevents the risk of water return to the device and minimizes condensation in the patient circuit to maintain a stable level of humidification overnight, for a therapy combining comfort and optimal efficiency.

#### **EASE OF USE**

The design features of DreamStar™ Intro help to facilitate the easy use of the device. They include automatic reminders, quiet operation and front screen illuminated retro buttons with automatic extinction of the backlight. The settings can be uploaded to the device via the DreamStar™ data card. Compliance data are stored for 365 days in the internal memory of the device for a simple analysis with DreamStar™ Analyze analysis software.

## **CPAP** device

#### **FEATURES:**

#### EASE OF USE:

- Front-facing screen
- · Auto reminder settings
- Large, backlit buttons
- · Multi language menu
- Can run on external battery (by-pass humidification option)
- Download compliance data via the DreamStar™ card
- Upload settings via the DreamStar<sup>™</sup> card
- Complete software analysis DreamStar™ Analyze

#### **TECHNOLOGIES:**

- DreamStar ™ Water Control
- Humid Control

#### **FULLY INTEGRATED:**

- Integrated HC heated humidifier
- Single electrical plug connection for both device and heated humidifier

#### **AUTOMATED:**

- Automatic altitude and barometric pressure change compensations
- Internal universal switching power supply
- · Automatic leak compensation

#### **SPECIFICATIONS:**

#### **Characteristics (or Performances):**

**PRESSURE OUTPUT**:  $4 \text{ cm H}_20 \text{ to } 20 \text{ cm H}_20 \text{ (+/- 0,6 cm H}_20)$ 

 ${\sf Max.}$  flow for  ${\sf 1cmH_2O}$  pressure drop at the patient output:

160 lpm at 20 cmH<sub>2</sub>O without humidifier

Noise LEVEL: 26 dB measured with an uncertainty of 2 dB according to

ISO 17510-1:2002

**OPERATIONAL MODE:** Constant pressure

RAMP DURATION: 0 to 45 minutes adjustable by 5 minute increments

COMFORT PRESSURE: 4 cm H<sub>2</sub>0 to prescribed pressure

**HUMIDIFIER CHAMBER CAPACITY:** 

300 ml at maximum water level

**HEATER ADJUSTABLE RANGE AND WATTAGE:** 

0 - 10 adjustable (40 Watts at setting 10)

**Physical Characteristics:** 

Size: 20 wide cm X 23 cm deep X 12 cm high

WEIGHT: 1,5 kg with humidifier; 1,4 kg without humidifier

FILTERS: Cleanable standard filter or disposable HEPA filter

**Power Characteristics:** 

Power SUPPLY: Input range 100-240 V, 50-60 Hz

BATTERY: 13 V

POWER CONSUMPTION:

50 VA max without humidifier; 90 VA with humidifier at setting 10

IEC 60601-1 classification : Class II (double isolation) Type BF

**Environmental Requirements:** 

**OPERATING TEMPERATURE:** 

5°C to 35°C, 10% to 95% relative humidity (RH), non-condensing

STORAGE TEMPERATURE:

-20°C to 60°C, relative humidity (RH) ≤95%, non-condensing

 $\textbf{Operating altitude}: \ from \ sea \ level \ 2700 \ meters$ 

#### **ORDERING INFORMATIONS:**

Products and accessories :	References:
DreamStar™ Intro (Includes patient circuit 1,8m and carrying case)	M-115500-UK
DreamStar™ Intro Evolve with reservoir (Includes patient circuit 1,8m and carrying case)	M-115520- UK
DreamStar™ Intro Evolve with cover (Includes patient circuit 1,8m and carrying case)	M-115510- UK
DreamStar™ Water Reservoir (includes base, cover and gasket)	M-315510-14
DreamStar ™Cover	M-315530-19
Filter holder kit (27x42)	M-314820-02
96 Inlet Air filter (27x42)	M-414840-06
100 High efficiency air filter (27x42)	M-414841-07
Carrying case DreamStar™	M-815505-00
Power cord 2 poles - UK	M-660601-40
Cigarette lighter cable 3A	M-214530-01
Battery cable 3A ring term	M-214530-00
Patient circuit 1,80 m / diam 22mm	M-261000-04
Patient manual DreamStar™ range	M-155DFU00-00
Practitioner manual DreamStar™ range	M-155DFU00-40

#### **Data Management:**

DreamStar™ Analyze software	M-215630-00
DreamStar™ Analyze clinical kit	M-215630-01
Memory card DreamStar™	M-215530-00
Memory card reader	M-214830-03
User manual	M-156DFU00-20



For more information contact us at the following addresses : customer service@sefam-medical.com marketing@sefam-medical.com

DreamStar™ is a trademark of SEFAM

€ 0459